





Filename: 2024-11-20 - 16.01.50.D5000

Sample Info

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
A1	72.9	Ladder		Caution! Expired ScreenTape device (used after two weeks of first use); Ladder
B1	22.1	PacuH10		Caution! Expired ScreenTape device (used after two weeks of first use)
C1	6.52	PcompE10		Caution! Expired ScreenTape device (used after two weeks of first use)


Compact Peak Table

Well	Sample Description	Size [bp]	Calibrated Conc. [ng/ul]	Assigned Conc. [ng/ul]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
A1	Ladder	15	7.96	-	816	-		Lower Marker
A1	Ladder	100	8.21	-	126	11.26		
A1	Ladder	250	9.21	-	56.7	12.62		
A1	Ladder	400	9.46	-	36.4	12.97		
A1	Ladder	600	10.2	-	26.1	13.97		
A1	Ladder	1000	10.3	-	15.9	14.19		
A1	Ladder	1500	9.43	-	9.67	12.92		
A1	Ladder	2500	7.33	-	4.51	10.05		
A1	Ladder	3500	5.42	-	2.38	7.43		
A1	Ladder	5000	3.36	-	1.03	4.60		
A1	Ladder	10000	3.25	3.25	0.500	-		Upper Marker
B1	PacuH10	15	6.62	-	679	-		Lower Marker
B1	PacuH10	234	22.1	-	145	100.00		
B1	PacuH10	10000	3.25	3.25	0.500	-		Upper Marker
C1	PcompE10	15	7.10	-	728	-		Lower Marker
C1	PcompE10	257	6.52	-	39.0	100.00		
C1	PcompE10	10000	3.25	3.25	0.500	-		Upper Marker



Default image (Contrast 100%)

A1: Ladder




Sample Table

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
A1	72.9	Ladder		Caution! Expired ScreenTape device (used after two weeks of first use); Ladder

Peak Table

Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	7.96	-	816	-		Lower Marker
100	8.21	-	126	11.26		
250	9.21	-	56.7	12.62		
400	9.46	-	36.4	12.97		
600	10.2	-	26.1	13.97		
1000	10.3	-	15.9	14.19		
1500	9.43	-	9.67	12.92		
2500	7.33	-	4.51	10.05		
3500	5.42	-	2.38	7.43		
5000	3.36	-	1.03	4.60		
10000	3.25	3.25	0.500	-		Upper Marker

B1: PacuH10




Sample Table

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
B1	22.1	PacuH10		Caution! Expired ScreenTape device (used after two weeks of first use)

Peak Table

Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	6.62	-	679	-		Lower Marker
234	22.1	-	145	100.00		
10000	3.25	3.25	0.500	-		Upper Marker

C1: PcompE10



Sample Table

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
C1	6.52	PcompE10		Caution! Expired ScreenTape device (used after two weeks of first use)

Peak Table

Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	7.10	-	728	-		Lower Marker
257	6.52	-	39.0	100.00		
10000	3.25	3.25	0.500	-		Upper Marker